

**Response To Pre-Bid Queries (Pre-Bid date: 24.04.2018)**

**NIB Ref: HITES/PCD/NCI-AIIMS/14/18-19 Dated: 13.04.2018**

**Schedule No. 24 (Ref/Event number 3000002800) VATS (Video Assisted Thoracoscopy Surgery) & Mediastinoscopy set with Full HD camera & Recording device**

Sl. No	Reference to the bidding document	Specification/Clause of the bidding documents	Name of firms	Queries/Clarification asked by firms	Recommendations/replies by the Committee
1	4.ii	Captures still HD images, video sequences (from 3 sources), and audio files	Karl Storz	Seek more clarification. We can record 2 sources simultaneously.	<b>Amended as:</b> Captures still HD images, video sequences (from 2 sources), and audio files
2	4.iii	Resolution of still images should be 1920x1080 and that of video 1080p @25p	Karl Storz	1080p 50 fps.	<b>Amended as:</b> Resolution of still images should be 1920x1080 and that of video 1080p
3	4.iv	Writes multi-session and multi-patient CDs/DVDs	Karl Storz	To be deleted	<b>Amended as:</b> Writes multi-session and multi-patient CDs/DVDs/USB
4	4.xii	HIPAA compliant	Karl Storz	To be deleted	Deleted
5	4.xv	xv. Chipset: Intel® 855GME + Intel® 6300ESB Embedded	Karl Storz	To be deleted	Deleted
6	4.xvi	xvi. Processor: Intel® Pentium® M 735	Karl Storz	To be deleted	Deleted
7	4.xvii	xvii. Graphics: Intel® Extreme Graphics 2 Controller onboard	Karl Storz	To be deleted	Deleted
8	4.xviii	xviii. Grabber-card: DVI-D, SDI, S-Video, Composite;	Karl Storz	Since THE platform is 2D, 3D & 4K; it has HDMI/DVI and DP interface so amend it accordingly.	Deleted
9	4.xix.	xix. Audio: AC97/DDS.1 onboard	Karl Storz	To be deleted	Deleted
10	4.xx	xx. Optional: AGP Slot onboard	Karl Storz	To be deleted	Deleted
11	4.xxvi.	I/O ports: 2 x PS/2, 2 x Serial, 3 x RJ45 (LAN), 4 x USB 2.0 (1 x Front), 3 x Audio (Line In, Line Out and Microphone), VGA;	Karl Storz	To be deleted	Deleted
12	4.xxviii.	D COM and HL7 interface - The DICOM 3 interface shall be installed to the system in order to allow the surgeon to view all the DICOM 3 images stored in the PACS system on a digital light box within the operating rooms. Furthermore, all intra operative images recorded can be sent via the DICOM 3 interface to the PACS system for further processing.	Karl Storz	please clarify whether you need to be inbuilt with recording system or software in separate PC. We recommend separate PC with software.	Deleted
13	4.1.j	CO2 bottle holder, max. diameter 155mm – 1no.	Karl Storz	May be deleted as CO2 bottle is not required..	Deleted
14		Quantity of VATS & Mediastinoscopy set with Full HD Camera & Recording Device: 2 nos. (as per tender requirement)			<b>Amended as:</b> Quantity of VATS & Mediastinoscopy set with Full HD Camera & Recording Device: 1 no.