

**Response To Pre-Bid Queries (Pre-Bid date: 24.04.2018)**

**NIB Ref: HITES/PCD/NCI-AIIMS/14/18-19 Dated: 13.04.2018**

**Schedule No. 15 - Rigid Cysto-Resectoscope (Rfx no. 3000002791)**

Sl. No	Reference to the bidding document	Specification/Clause of the bidding documents	Name of firms	Queries/Clarification asked by firms	Recommendations/replies by the Committee
1	Page 82 Para 2a	Telescope of 4mm with 0, 30, 70 deg direction of view			Telescope of 4mm with 0, 30, 70 deg direction of view - 1 each
2	Page 83 BOQ Sl.no.1	Telescope of 4mm with 0, 30, 70 deg direction of view - 1 no.			Telescope of 4mm with 0, 30, 70 deg direction of view - 1 no. each
3	Page 82 Para 3	300 Watt Xenon/LED Light Source	Olympus	Kindly add Light guide cable with length of 4 mm or more	Added Para: Light guide cable 2 nos. with length of 3mtr or more should be provided along with light source
4			Stryker	Endovision system not asked in the tender	No Change
5		Kindly add Bipolar Cautery as same has not been added for Resection	Olympus	Kindly add technical specs for Bipolar cautery: Electro-Surgical Unit: 300 Watt or more  <b>a. Microprocessor Controlled Electrosurgical Unit: Should have following specifications:</b> i. Should have built-in regulated fan ii. Should have LCD and Touch Screen user Interface iii. Should have memory function with characters iv. Should have Automatic Saline Detection mode  <b>b. Monopolar Cut Modes :</b> i. PureCut for continuous cutting ii. BlendCut for Continuous cutting with Coagulation iii. Cut and coagulation modes for Intermittent cutting with long coagulation  <b>c. Monopolar Coagulation Modes :</b> i. SoftCoag or similar mode for Soft and deep coagulation ii. PowerCoag or similar mode for Fast coagulation and dissection iii. ForcedCoag or similar mode for Fast coagulation iv. SprayCoag or similar mode for Contact-free superficial coagulation  <b>d. Bipolar Cut Modes :</b> i. BipolarCut or similar mode for Continuous cutting ii. SalineCut or similar mode for Continuous cutting in saline  <b>e. Bipolar Coagulation Modes:</b> i. BiSoftCoag or similar mode for Soft coagulation (+/- Autostart) ii. AutoCoag or similar mode for Soft coagulation with Autostop iii. SalineCoag or similar mode for Coagulation in saline iv. FineCoag or similar mode Fine coagulation	No Change

*[Handwritten signatures and initials in blue ink]*



Sl. No	Reference to the bidding document	Specification/Clause of the bidding documents	Name of firms	Queries/Clarification asked by firms	Recommendations/replies by the Committee
				v. RFCoag for or similar mode Deep (tumor) ablation with Autostop	
6	Page 82 Para 1.c	Should have Robust ceramic beak with lifetime guarantee against thermal damage	Karl Storz	Kindly remove word "life time guarantee" for ceramic beak.	Amended as: Should have Robust ceramic beak with guarantee during warranty period
7	Page 82 Para 2.n	Compatible Monopolar and Bipolar HF Cable, 4 m to be provided	Karl Storz	Kindly amend the length in a range of 3-4 m.	Amended as: Compatible Monopolar and Bipolar HF Cable, 3 mtr or more to be provided
8	Page 83 BOQ Sl.no.14	Compatible Monopolar and Bipolar HF Cable, 4 m to be provided - 1no.			Amended as: Compatible Monopolar and Bipolar HF Cable, 3 mtr. or more to be provided - 1 no.
9	Page 82 Para 2.t	Bipolar HF-Resection Electrode, Roller type or similar (5 pc)	Karl Storz	Kindly change this to Roller type/Half moon shape, 6pcs.	Amended as: Bipolar HF-Resection Electrode, Roller type or halfmoon shape or similar (5 pc)
10	Page 83 BOQ Sl.no.20	Bipolar HF-Resection Electrode, Roller type or similar -5 nos.			Amended as: Bipolar HF-Resection Electrode, Roller type or halfmoon shape or similar - 5 nos.
11	Page 82 Para 2.u	BipolarHF-Resection Electrode, 45 Deg, Needle or similar (6 pc)	Karl Storz	Kindly remove word 45 degree and it should be Pointed coagulation electrode.	Amended as: BipolarHF-Resection Electrode, Needle or similar (6 pc)
12	Page 83 BOQ Sl.no.21	BipolarHF-Resection Electrode, 45 Deg, Needle or similar - 6 nos.			Amended as: BipolarHF-Resection Electrode, Needle or similar - 6 nos.
13	Page 82 Para 2.v	BipolarHF-Resection Electrode, Medium Cutting Loop or similar (12 pc)	Karl Storz	Kindly amend this to medium or 24/26 cutting loop bipolar.	Amended as: BipolarHF-Resection Electrode, medium or 24/26 Cutting Loop or similar (12 pc)
14	Page 83 BOQ Sl.no.22	BipolarHF-Resection Electrode, Medium Cutting Loop or similar - 12 nos.			Amended as: BipolarHF-Resection Electrode, medium or 24/26 Cutting Loop or similar - 12 nos.
15	Page 82 Para 2.x	BipolarHF-Resection Electrode, Button, for plasma vaporization or similar (12 pc)	Karl Storz	Kindly amend this to plasma vaporization/vapo-enucleation process.	Amended as: BipolarHF-Resection Electrode, Button, for plasma vaporization/vapo-enucleation process or similar (12 pc)
16	Page 83 BOQ Sl.no.24	BipolarHF-Resection Electrode, Button, for plasma vaporization or similar -12 nos.			Amended as: BipolarHF-Resection Electrode, Button, for plasma vaporization/vapo-enucleation process or similar - 12 nos.
17	Page 82 Para 2.y	Monopolar HF-Resection Electrode, cutting (5 pcs)	Karl Storz	Kindly amend this to Angled cutting loop, 6pcs.	Amended as: Monopolar HF-Resection Electrode, cutting /angled cutting loop - 5 nos.



Sl. No	Reference to the bidding document	Specification/Clause of the bidding documents	Name of firms	Queries/Clarification asked by firms	Recommendations/replies by the Committee
	Page 83 BOQ Sl.no.25	Monopolar H-F-Resection Electrode, cutting - 5 nos.			<b>Amended as:</b> Monopolar HF-Resection Electrode, cutting /angled cutting loop - 5 nos.
18	Page 82 Para 2.z	Monopolar HF-Resection Electrode, coagulation (5 pcs)	Karl Storz	Kindly amend this to Coagulation electrode/ball shaped, 6 pcs	<b>Amended as:</b> Monopolar HF-Resection Electrode, coagulation/ball shaped - 5 nos.
19	Page 83 BOQ Sl.no.26	Monopolar HF-Resection Electrode, coagulation - 5 nos.			<b>Amended as:</b> Monopolar HF-Resection Electrode, coagulation/ball shaped - 5 nos.
20	Page 82 Para 3	300 Watt Xenon/LED Light Source :	Karl Storz	Light Source is already offered at many urological items: Do you really need this light source again? If yes then below changes are required.	<b>No Change</b>
21	Page 82 Para 3.a	Separate/combined module from video processor	Karl Storz	Delete	<b>Deleted</b>
22	Page 82 Para 3.c	Powerful 300 Watt Xenon Lamp with backup lamp of Halogen/LED	Karl Storz	The back up lamp of Halogen/LED is a very company specific and usually it is not required with powerful xenon 300 watt light source therefore kindly delete back up lamp of Halogen/LED.	<b>Amended as:</b> Powerful 300 Watt Xenon Lamp/LED light source.
	Page 82 Para 3.d	Lamp life should not be less than 500 hrs (approx), one spare 300W xenon lamp should be supplied in case of xenon light source			<b>Amended as:</b> Lamp life should not be less than 500 hrs (approx). 2 nos. spare 300W xenon bulbs should be supplied in case of xenon light source
					<b>Added Para under light source:</b> Should be supplied with Light guide cables 3 mtr or more length - 2 nos. with each light source.  Adaptors for connecting to rigid scopes (proctoscope, sigmoidoscope etc.) of other reputed makes should be offered
23	Page 82 Para 4	Supplier should provide OEM Trolley along with the system	Karl Storz	Delete	<b>Amended as:</b> Supplier should provide Trolley along with the system from same principle manufacturer

*Sanjay*  
*Dr. Sanjay*

*Dr.*  
*Dr.*

